

HI-CREPE COTTON CREPE BANDAGE

Features of the Product:

- Superior flesh coloured elastic bandage with non-fraying selvedge gives uniform constant elasticity, combination of cotton stretch material with excellent dermophilic properties.

Indications:

Post treatment for all fracture dislocations,
Sprains and strains, Skin grafting

Sizes:

5 cm x 4 m
8 cm x 4 m
10 cm x 4 m
15 cm x 4 m
20 cm x 7 m



References:

Melhuish, J., et al., The Physics of sub-bandage pressure measurement. Journal of Wound Care, 2000. 9(7): p. 308-310.

Fentem PH, Goddard M, Gooden BA. Support for varicose veins. Br Med J. 1976 Jan 31; 1(6004):254-256.

[PMC free article] [PubMed]

Callam, M.J.Harper, D.R. Dale, J.J et al, Lothian and Forth Valley leg ulcer healing trial.

Part 1. Elastic versus non-elastic bandaging in the treatment of chronic leg ulceration. Phlebology. 1992; 7:130

Eriksson, G. Comparison of two occlusive bandages in the treatment of venous leg ulcers. Br J Dermatol. 1986;114:227-230

Lippmann, M.D. Fishman, L.M,Farrar, R.H et al, Edema control in the management of disabling chronic venous insufficiency. Arch Phys Med Rehabil. 1991;75:436-441

ABOUT US:

We are a medical device manufacturing and associated service providing company. Established with a vision to provide a world class medical experience, second to none, we at OUTCARE are guided by hard-core business principles hand in glove with our commitment for the betterment of healthcare sector.



Manufactured & Marketed by :
OUTCARE MEDICAL DEVICES PVT LTD
Second Floor, Amit Plaza, Somajiguda Main Circle,
Hyderabad, Telangana, INDIA - 500 082

For Technical Information write us at info@outcaremed.com
Customer Care : +91 93989 24791
Visit us on www.outcaremed.com



COMPRESSION BANDAGES

HI-GRIP VC ADVANCED COMPRESSION BANDAGE

Features of the Product:

- Advance Compression bandage with a pocket for easy fit and quick wrapping.
- Adjustable Velcro strap for comfort and structural support.
- No clips or tapes required to hold the bandage in place.
- Premium quality breathable material with 200% extensibility.
- Excellent regain ability ensures consistent and controlled compression.
- High quality material retains longevity and compression even after multiple washes.

Sizes:

10 cm × 4 m
15 cm × 4 m
15 cm × 8 m



HI-GRIP VC COMPRESSION BANDAGE PRESSURE STANDARDS

The compression pressure was measured by half overlap and the average readings are given below:

HALF OVERLAP					
Position of measurement	Above ankle (b)	Above ankle position (b1)	Calf Position (c)	Below knee position (d)	Above knee position (f)
Pressure unit	Average mmHg	Average mmHg	Average mmHg	Average mmHg	Average mmHg
Elastic Compression Bandage	28.86	27.4	25.8	22.8	18.6

The bandage was tested for compression pressure by half overlap with uniform stretch from toe to knee. The pressure of 25-29mmHg more consistently when applying single layer bandage than with the double layer bandage.

Indications:

- Venous Disorders – Chronic Venous Insufficiency, Venous, Ulcers, Varicose Veins, Deep Vein Thrombosis, Lymphatic Disorders, Post-operative edema – CABG, Post Total Knee Replacement

WRAPPING INSTRUCTIONS VIDEO <https://www.youtube.com/watch?v=6Gs7xcnCDXo&t=12s>



HI-GRIP PREMIUM (Cotton Compression Bandage)

Features of the Product:

- Superior Stretching with more than 200% extensibility.
- Excellent regain ability ensures consistent and controlled compression.
- Conforms to body contour.
- Hypoallergenic and skin friendly for sensitive skin.

Sizes:

10 cm × 4 m
15 cm × 4 m
15 cm × 8 m

Indications:

Lymphatic disorders, Venous Disorders, Post Phelbectomy, Lymphoedema



FOURBAN 4 LAYER COMPRESSION BANDAGE

Description

Fourban - 4 layer compression bandage kit consists of 4 compression wraps.

Includes

Layer 1 – Orthopaedic bandage provides cushioning effect and protects the bony prominence.
 Layer 2 – Conforming bandage protects the bony prominence with mild compression.
 Layer 3 – Compression bandage provides the graduated compression from ankle to the calf and minimizes the edema and increases venous return. (Figure 8 application)
 Layer 4 – Cohesive bandage provides effective and sustained compression throughout the wear time.

Features of the Product

- Graduated compression of lower leg aids venous blood circulation.
- Absorbent for exudate management.
- Fewer dressing changes, wear time of dressing is upto seven days.
- Cost effective.

Layer	Type of Bandage	Size
Layer 1	Orthopaedic Bandage	10cm x 4m (Unstretched)
Layer 2	Conforming Bandage	10cm x 4.5m (Stretched)
Layer 3	Compression Bandage	10cm x 9.14m (Stretched)
Layer 4	Cohesive Bandage	10cm x 4.5m (Stretched)

Indications:

CVI-Chronic Venous Insufficiency (with/without wound), Lymphatic Disorders (Edema, edematous limbs), Elephantiasis



WRAPPING INSTRUCTIONS VIDEO <https://www.youtube.com/watch?v=z3Y4mmayGdo>



HI-GRIP ELASTIC COMPRESSION BANDAGE

Features of the Product:

- Superior Stretching properties.
- Excellent regain ability ensures consistent and controlled compression.
- More than 200% extensibility.
- The pressure is between 25-29mm of Hg.

Sizes:

8 cm × 4 m
10 cm × 4 m
15 cm × 4 m

Indications:

Post treatment for all fracture dislocations, Sprains and strains, Skin grafting, Varicose veins

