

### Who We Are



Bio Compression has been the leading manufacturer of pneumatic compression devices for almost 40 years.

All of our pumps are lightweight and user friendly with a non-peristaltic action that mimics the body's own natural function. This action helps to increase circulation and reduce edema to manage your symptoms and improve your quality of life.

Compression therapy is used for the treatment of lymphedema, lipedema, venous insufficiency, wound healing, DVT prophylaxis, and many other circulatory disorders.

Based in Moonachie, New Jersey – our manufacturing facility is focused on quality products, boasting the largest selection of garments and the longest guarantee in the industry.

With countless clinical studies to prove the efficacy of our products, it's no wonder that patients from all over the world use Bio Compression devices to help treat their lymphatic and circulatory disorders. We are proud to help you.













### 2000 SERIES

### On the forefront of pneumatic compression is our 2000 series pumps. Designed to be user friendly, a digital

screen allows for easy operation. The pause button allows you to stop treatment and start again without skipping a beat; focus therapy allows for extra attention to a specific area; and the compliance meter keeps track of your hours of use. These gradient sequential pneumatic compression devices are used for the treatment of lymphedema, venous insufficiency, and various other circulatory disorders. Improve your quality of life starting today with the 2000 series pumps.





### SC-2008-DL

### Benefits Include:

- Manage pain
- Increase lymphatic and vascular flow
- Promote wound healing
- Reduce swelling
- Prevent infection
- Improve quality of life

### **Indications**

- Lymphedema
- Lipedema
- Phlebolymphedema
- DVT prophylaxis
- · Chronic venous insufficiency
- Wound care

### Key Features:

SC-2004-DL

- 4-Chamber compression therapy
- Unilateral or bilateral capability
- Adjustable pressure between 10 mmHg and 120 mmHg
- Adjustable cycle time between 60 and 120 seconds
- Countdown timer or continuous mode
- Focus therapy
- Pause button
- Compliance meter
- 3 year warranty

Specifications:

Dimensions: 4.5" x 12" x 7.34"

11 cm x 30 cm x 18 cm

Weight: 3.3 lbs. 1,5 kgs

HCPCS: E0651

### **Key Features:**

- 8-Chamber compression therapy
- Unilateral or bilateral capability Adjustable pressure between
- Adjustable pressure between 10 mmHg and 120 mmHg
- Adjustable cycle time between 60 and 120 seconds
- Countdown timer or continuous mode
- Focus therapy
- Pause button
- Compliance meter
- 3 year warranty

Specifications:

Dimensions: 4.5" x 12" x 7.34"

11 cm x 30 cm x 18 cm

Weight: 3.65 lbs. 1,7 kg

HCPCS: E0651





See product manual for warnings, precautions, and contraindications.


# Made to the hig



Bio Compression Systems produces the most extensive range of garments in the industry. With over 50 different standard sizes available off the shelf for lower extremity, upper extremity, abdominal and truncal patients, nobody is left untreated by Bio Compression. And for those rare occasions where our standard sleeves do not fit, our manufacturing facility produces custom garments once a week to assure everyone is properly fitted without having to wait.

All of our garments are made of high quality 200 Denier Nylon Oxford, coated with 3 mil of polyurethane to assure a level of water resistance that allows all of our sleeves to be cleaned. Engineered with our special baffle system to guarrantee there is no gapping during treatment, our sleeves are latex free and do not require expanders, separate tubing, or any other accessories. A simple connection system allows you to plug the sleeves into our pump with ease and be operational within seconds.

All of our lower extremity sleeves come complete with pull-up straps for ease of use, and our entire range of sleeves includes protective loops to stop unnecessary pulling of the tubes. Biocomfort Garments are manufactured with zippers from YKK, the best of its kind. We spare no cost in assuring that you get the greatest sleeve possible. Our goal is 100% satisfaction for every patient we treat.



### Baffle System

While many other companies promote an overlapping chamber system with their sleeves, our garments have been carefully engineered with a baffle system that assures the best flow of pressure from one chamber to the next. The construction of independently controlled pressure chambers assures there is never gapping, peaks and valleys or bleed-through with the delivered pressure.

### hest Standards




### Quick Connect System

Our Biocomfort Garments all come ready to use. Just line up the chamber numbers on our specially developed clips with the numbers on our pumps, squeeze the clip and insert. It's that easy to get your treatment sessions started.



### YKK<sup>®</sup> Zippers

We spare no cost in assuring that your garment is the best possible product. All of our zippers come from YKK®, one of the leaders in the industry. Complete with automatic locking sliders, our sleeves will not unzip accidently when inflated.



### Welded Air Connection Ports

All of our sleeves come complete with connection ports that are RF Welded into the sleeves. This construction guarantees that there will be no leaking, no breaking, and a perfect flow of air each and every time you sit down for a treatment.



### Integrated Pull-Up Straps


All of our lower extremity garments come complete with built-in pull-up straps. While our garments are already the easiest to don and doff, we have taken every step possible to assure your patients are compliant every step of the way.

## 4 Chambe STANDARD LEG SLEEVE





	4 Chamber Standard Leg Sleeve Size and Part Number										
		Measur	ements belo	w are in <b>Inch</b>	ies and are th	e maximum	allotted poin	t on limb			
				Circui	mferences				Lengths		
Part Number	Si	ize	Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease	
GN-3045-S	Small	Narrow	25.5"	25"	22*	20"	18"	19.5"	25.5"	27"	
GN-3045-M	Medium	Narrow	27.5"	25.5"	23*	20"	18"	19.5"	28"	32"	
GN-3045-L	Large	Narrow	28"	26"	24*	21"	18"	19.5"	32"	36"	
GS-3045-H	Half	Standard	N/A	N/A	29"	27"	23.5"	22"	19"	NA	
GS-3045-S	Small	Standard	31"	30"	28"	24.5"	21.5"	19.5"	25.5"	27"	
GS-3045-M	Medium	Standard	32"	30"	28"	24.5"	21.5"	19.5"	28"	32"	
GS-3045-L	Large	Standard	33"	32"	30"	26"	21.5"	19.5"	32"	36"	
GW-3045-H	Half	Wide	N/A	N/A	30.5"	28"	29"	23"	19"	NA	
GW-3045-S	Small	Wide	35.5"	34"	30.5*	28"	25.5"	21.5"	25.5*	27"	
GW-3045-M	Medium	Wide	36.5"	34"	30"	28"	25.5"	22.5"	28"	32"	
GW-3045-L	Large	Wide	37*	36"	32"	29"	24"	20.5"	32"	36"	
GXW-3045	Medium	Extra Wide	44"	41"	36"	34"	30"	22.5"	29"	32"	
		Measurem	ents below a	re in <b>Centim</b>	eters and are	the maximu	ım allotted p	oint on limb			
				Circui	mferences				Lengths		
GN-3045-S	Small	Narrow	65	64	56	51	46	50	65	69	
GN-3045-M	Medium	Narrow	70	65	58	51	46	50	71	81	
GN-3045-L	Large	Narrow	71	66	61	53	46	50	81	91	
GS-3045-H	Half	Standard	N/A	N/A	74	69	60	56	48	NA	
GS-3045-S	Small	Standard	79	76	71	62	55	50	65	69	
GS-3045-M	Medium	Standard	81	76	71	62	55	50	71	81	
GS-3045-L	Large	Standard	84	81	76	66	55	50	81	91	
GW-3045-H	Half	Wide	N/A	N/A	78	71	74	58	48	NA	
GW-3045-S	Small	Wide	90	86	77	71	65	55	65	69	
GW-3045-M	Medium	Wide	93	86	76	71	65	57	71	81	
GW-3045-L	Large	Wide	94	91	81	74	61	52	81	91	
GXW-3045	Medium	Extra Wide	112	104	91	86	76	57	74	81	



HCPCS: E0669 - Half leg E0667 - Full leg

## 4 Chambe ADJUSTABLE LEG SLEEVE






### 4 Chamber Adjustable Leg Sleeve Size and Part Number

				Circur	mferences				Lengths		
Part Number	Size		Thigh Mid Mid Mid Ankle Arch below Thigh Knee Calf (Instep Crotch		Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease				
GA-3045-H	Half	Adjustable	N/A	N/A	32"	30"	23.5"	22"	19*	NA	
GA-3045-S	Small	Adjustable	34"	33"	31"	27.5"	21.5"	19.5"	25.5"	27"	
GA-3045-M	Medium	Adjustable	35"	34"	31"	27.5"	21.5"	19.5"	28*	32"	
GA-3045-L	Large	Adjustable	37"	36"	33"	29"	21.5"	19.5"	32"	36"	
GWA-3045-H	Half	Adj. Wide	N/A	N/A	33.5"	31"	29"	23"	19"	N/A	
GWA-3045-S	Small	Adj. Wide	38.5"	37"	33.5"	31"	25.5"	21.5"	25.5"	27"	
GWA-3045-M	Medium	Adj. Wide	39.5"	36.5"	33.5"	31"	25.5"	22.5"	28"	32"	
GWA-3045-L	Large	Adj. Wide	40"	38"	35"	32"	24"	22.5"	32"	36"	
GXWA-3045	Medium	Adj. Extra Wide	47"	44"	39"	37"	30"	22.5"	29"	32"	

Measurements below are in (	Centimeters and are the maximu	m allotted point on limb
Treasurements setom are mig	and are the manning	m accorded points on time

				Circu	mferences				Lengths	
GA-3045-H	Half	Adjustable	N/A	N/A	81	76	60	56	48	NA
GA-3045-S	Small	Adjustable	86	84	79	70	55	49	65	69
GA-3045-M	Medium	Adjustable	89	86	79	70	55	49	71	81
GA-3045-L	Large	Adjustable	94	91	84	74	55	49	81	91
GWA-3045-H	Half	Adj. Wide	N/A	N/A	85	79	74	58	48	NA
GWA-3045-S	Small	Adj. Wide	98	94	85	79	65	55	65	69
GWA-3045-M	Medium	Adj. Wide	100	93	85	79	65	57	71	81
GWA-3045-L	Large	Adj. Wide	102	96	89	81	61	57	81	91
GXWA-3045	Medium	Adj. Extra Wide	119	112	99	94	76	57	74	81



HCPCS: E0669 – Half leg E0667 – Full leg



## 8 Chambe STANDARD LEG SLEEVE






### 8 Chamber Standard Leg Sleeve Size and Part Number Measurements below are in Inches and are the maximum allotted point on limb

				Circu	mferences				Lengths			
Part Number	Size		Size		Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease
G8-3045-S	Small	Standard	31"	30"	28"	24.5"	21.5"	19.5"	25.5*	27"		
G8-3045-M	Medium	Standard	32"	30"	28"	24.5"	21.5"	19.5"	28"	32"		
G8-3045-L	Large	Standard	33"	32"	30"	26"	21.5"	19.5"	32"	36"		
GN8-3045-S	Small	Narrow	25.5"	25"	22"	20"	18"	19.5"	25.5"	27"		
GN8-3045-M	Medium	Narrow	27.5"	25.5"	23"	20"	18"	19.5"	28"	32"		
GN8-3045-L	Large	Narrow	28"	26"	24"	21"	18"	19.5"	32"	36"		
GW8-3045-S	Small	Wide	35.5"	34"	30.5"	28"	25.5"	21.5"	25.5*	27"		
GW8-3045-M	Medium	Wide	36.5"	34"	30"	28"	25.5"	22.5"	28"	32"		
GW8-3045-L	Large	Wide	37"	36"	32"	29"	24"	20.5"	32"	36"		
GXW8-3045	Medium	Extra Wide	44"	41"	36"	34"	30"	22.5"	29"	32"		

Measurements below are in <b>Centimeters</b> and are the maximum allotted point on limb
---

				Circui	nferences				Lengths				
G8-3045-S	Small	Standard	79	76	71	62	55	49	65	69			
G8-3045-M	Medium	Standard	81	76	71	62	55	49	71	81			
G8-3045-L	Large	Standard	84	81	76	66	55	49	81	91			
GN8-3045-S	Small	Narrow	65	63	56	51	46	49	65	69			
GN8-3045-M	Medium	Narrow	70	65	58	51	46	49	71	81			
GN8-3045-L	Large	Narrow	71	66	61	53	46	49	81	91			
GW8-3045-S	Small	Wide	90	86	77	71	65	55	65	69			
GW8-3045-M	Medium	Wide	93	86	76	71	65	57	71	81			
GW8-3045-L	Large	Wide	94	91	81	74	61	52	81	91			
GXW8-3045	Medium	Extra Wide	112	104	91	86	76	57	74	81			

HCPCS: E0667 - Full leg




## 8 Chamber ADJUSTABLE LEG SLEEVE





8	8 Chamber Adjustable Leg Sleeve Size and Part Number										
		Measur	ements belov	v are in <b>Inch</b>	es and are th	e maximum a	allotted point	on limb			
				Circu	mferences				Lengths		
Part Number	Si	ze	Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease	
A8-3045-S	Small	tall Adjustable 34" 33" 31" 27.5" 21.5" 19.5"						19.5"	25.5"	27"	
A8-3045-M	Medium	Adjustable	35"	34"	31"	27.5"	21.5"	19.5"	28*	32"	
A8-3045-L	Large	Adjustable	37"	36"	33"	29"	21.5"	19.5"	32"	36"	
GWA8-3045-S	Small	Adj. Wide	38.5"	37"	33.5"	31"	25.5"	21.5"	25.5"	27"	
GWA8-3045-M	Medium	Adj. Wide	39.5"	36.5"	33.5"	31"	25.5"	22.5"	28"	32"	
GWA8-3045-L	Large	Adj. Wide	40"	38"	35"	32"	24"	22.5"	32"	36"	
GXWA8-3045	Medium	Medium Adj. Extra Wide 47" 44" 39" 37" 30" 22.5"		22.5"	29*	32"					
		Measurem	ents below a	e in <u>Centim</u>	eters and are	the maximu	n allotted po	int on limb			
				Circu	mferences				Lengths		
A8-3045-S	Small	Adjustable	86	84	79	70	55	49	65	69	
A8-3045-M	Medium	Adjustable	89	86	79	70	55	49	71	81	
A8-3045-L	Large	Adjustable	94	91	84	74	55	49	81	91	
GWA8-3045-S	Small	Adj. Wide	98	94	85	79	65	55	65	69	
GWA8-3045-M	Medium	Adj. Wide	100	93	85	79	65	57	71	81	
GWA8-3045-L	Large	Adj. Wide	102	96	89	81	61	57	81	91	
GXWA8-3045	Medium	Adj. Extra Wide	119	112	99	94	76	57	74	81	

HCPCS: E0667 - Full leg




### 4 Chamber





4 (	4 Chamber Arm Sleeve Size and Part Number										
	Measurements below are in <u>Inches</u> and are the maximum allotted point on limb										
Circumferences Lengths											
Part Number	Size	Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips			
GS-3035-S	Small	23.5"	22"	19.5"	17.5"	15"	24"	26.5"			
GS-3035-M	Medium	27"	26"	24"	21"	18"	28"	31"			
GS-3035-L	Large	25.5"	24"	21.5"	19"	16"	31"	35.5"			
M	1easurement	s below are ir	Centimeter	<b>s</b> and are the	maximum al	lotted point	on limb				
			Cir	rcumferences	;		Len	gths			
GS-3035-S	Small	60	56	49	44	38,	61	67			
GS-3035-M	Medium	69	66	61	53	46	71	79			
GS-3035-L	Large	65	61	55	48	41	79	90			

HCPCS: E0668



### 4 Chamber Shoulder Sleeve





### 4 Chamber Arm & Shoulder Sleeve Size and Part Number Measurements below are in $\underline{\textbf{Inches}}$ and are the maximum allotted point on limb Circumferences Lengths Axilla Bicep Elbow Forearm Wrist Axilla to Shoulder to Part Number Size Fingertips Fingertips GS-3035-SH-S-L Small Left 23.5 22' 19.5" 17.5" 15" 24" 26.5"+8" GS-3035-SH-S-R Small Right 23.5" 22" 19.5" 17.5" 15" 24" 26.5"+8" GS-3035-SH-M-L 27" 31"+8" Medium Left 26" 24" 21" 18" 28" GS-3035-SH-M-R Medium Right 27" 26" 24" 21" 18" 28" 31"+8" GS-3035-SH-L-L 25.5 24" 21.5" 19" 16" 31" 35.5"+8" Large Left GS-3035-SH-L-R Large Right 25.5" 24" 21.5" 19" 16" 31" 35.5"+8" Measurements below are in **Centimeters** and are the maximum allotted point on limb Circumferences Lengths GS-3035-SH-S-L Small Left 56 38 61 67+20 GS-3035-SH-S-R Small Right 67+20 60 56 49 44 38 61 GS-3035-SH-M-L Medium Left 66 61 53 46 71 79+20 GS-3035-SH-M-R 61 53 46 71 79+20 Medium Right 69 66

61

55

48

41


79

90+20 90+20

GS-3035-SH-L-R HCPCS: E0668

Large Right

GS-3035-SH-L-L




## 8 Chamber





8 (	8 Chamber Arm Sleeve Size and Part Number										
	Measurements below are in <u>Inches</u> and are the maximum allotted point on limb										
	Circumferences Lengths										
Part Number   Size   Tantan   Tantan								Shoulder to Fingertips			
G8-3035-S	Small	23.5"	22"	19.5"	17.5"	15*	24"	26.5"			
G8-3035-M	Medium	27"	26"	24"	21"	18"	28"	31"			
G8-3035-L	Large	25.5"	24"	21.5"	19"	16"	31"	35.5"			
	<u> </u>										
Ŋ	1easurement	s below are ir	<u>Centimeter</u>	<b>s</b> and are the	maximum al	lotted point	on limb				
			Cir	rcumferences	i		Len	gths			
G8-3035-S	Small	60	56	49	44	38	61	67			
G8-3035-M	Medium	69	66	61	53	46	71	79			
G8-3035-L	Large	65	61	55	48	41	79	90			

HCPCS: E0668



### ARM & SHOULDER SLEEVE





### 8 Chamber Arm & Shoulder Sleeve Size and Part Number Measurements below are in <u>Inches</u> and are the maximum allotted point on limb Circumferences Lengths

Part Number	Size	Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips
G8-3035-SH-S-L	Small Left	23.5"	22"	19.5*	17.5"	15"	24"	26.5"+8"
G8-3035-SH-S-R	Small Right	23.5"	22"	19.5*	17.5"	15"	24"	26.5"+8"
G8-3035-SH-M-L	Medium Left	27"	26"	24"	21"	18"	28"	31"+8"
G8-3035-SH-M-R	Medium Right	27"	26"	24"	21"	18"	28"	31"+8"
G8-3035-SH-L-L	Large Left	25.5"	24"	21.5*	19"	16"	31"	35.5"+8"
G8-3035-SH-L-R	Large Right	25.5"	24"	21.5*	19"	16"	31"	35.5"+8"
	•	•	•	•		•		

	Measurements	below are ir	<u>Centimeter</u>	<b>s</b> and are the	maximum al	lotted point	on limb			
Circumferences Lengths										
G8-3035-SH-S-L	Small Left	60	56	49	44	38	61	67+20		
G8-3035-SH-S-R	Small Right	60	56	49	44	38	61	67+20		
G8-3035-SH-M-L	Medium Left	69	66	61	53	46	71	79+20		
G8-3035-SH-M-R	Medium Right	69	66	61	53	46	71	79+20		
G8-3035-SH-L-L	Large Left	65	61	55	48	41	79	90+20		
G8-3035-SH-L-R	Large Right	65	61	55	48	41	79	90+20		

HCPCS: E0668

