

# co-wraps™

## INELASTIC COMPRESSION WRAP




## Say Goodbye to Swelling!

Our inelastic compression wrap improves venous and lymphatic return, effectively reducing edema and helping you feel better, faster.

### Features of the Product

- Inelastic compression wrap with pocket for easy fit and quick wrap.
- Perforated, breathable material that is gentle on the skin.
- The wrap consists of series of inelastic bands with velcro strap for comfort and structural support.
- Easy for donning and doffing with regular and criss cross wrap pattern.
- Inelastic wrap with adjustable pressure of 20-50mmHg.




### SIZE CHART (B-D)


L	Small	Medium	Large	XLarge	XXLarge
D	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
C	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
B1	24-34cm	29-39cm	34-44cm	39-49cm	44-54cm
B	19-26cm	21-30cm	26-36cm	31-41cm	36-46cm

## CALF


WRAP. RELIEVE. RESTORE.




Pocket for easy fit and quick wrap



Regular and Criss Cross Wrap Pattern



Adjustable pressure of 20-50mmHg



[WWW.OUTCAREMED.COM](http://WWW.OUTCAREMED.COM)

### FOR TECHNICAL INFORMATION

WRITE US AT [INFO@OUTCAREMED.COM](mailto:INFO@OUTCAREMED.COM)  
VISIT US ON [WWW.OUTCAREMED.COM](http://WWW.OUTCAREMED.COM)  
CUSTOMER CARE: +9140 48593542, +91 9398924791  
NOT MADE WITH NATURAL RUBBER LATEX



# co-wraps™

## INELASTIC COMPRESSION WRAP




## Your Recovery Starts Here.

Our inelastic compression wrap improves venous and lymphatic return, effectively reducing edema and helping you feel better, faster.

### Features of the Product

- Inelastic compression wrap with comfort and structural support.
- Perforated, breathable material that is gentle on the skin.
- The wrap consists of series of inelastic bands and strap extender with velcro strap to secure it firmly.
- Easy for donning and doffing with regular and criss cross wrap pattern.
- Inelastic wrap with adjustable pressure of 20-50mmHg.

### DETERMINING THE CORRECT SIZE




- C** the distance from ankle joint to the tip of the toe.
- A** The circumference of foot measured around its widest part.
- B** The circumference above foot measured around its narrowest part.
- Compare the measured values against the table below

### SIZE CHART


	XSmall	Small	Medium / Large
A	22-28cm	28-36cm	≤39cm
B	18-22cm	22-26cm	≤30cm
C	10-14cm	14-18cm	18-20cm

## FOOT


WRAP. RELIEVE. RESTORE.




Breathable material that is gentle on skin



Strap Extender with velcro strap



Adjustable pressure of 20-50mmHg



WWW.OUTCAREMED.COM

### FOR TECHNICAL INFORMATION

WRITE US AT [INFO@OUTCAREMED.COM](mailto:INFO@OUTCAREMED.COM)  
VISIT US ON [WWW.OUTCAREMED.COM](http://WWW.OUTCAREMED.COM)  
CUSTOMER CARE: +9140 48593542, +91 9398924791  
NOT MADE WITH NATURAL RUBBER LATEX

