

Comfortable, Easy to wear Cotton Stockings

The Classic Line of cotton Stockings, Particularly suitable for patients who want a Soft contact with skin & an excellent Quality Vs Price ratio.

Thanks for the structure of the Stockings, it ensures perfect transpiration And minimum possibility of skin allergies. Ideal for over-the-counter sales.







COMPRESSION CLASSES:

VEGAL 1: 15-20 mmHg VEGAL 2: 22-28 mmHg

INDICATIONS

15-20 mmHg in the presence of light venous disorders such as heaviness and a sense of fatigue in the legs, in non-serious post traumatic situations, after removal of cast, to support the muscles.

22-28 mmHg light CVI and after: varicophlebitis, healing of a small superficial ulcer, sclerosing treatment, selective outpatient variceal surgery, stripping, more pronounced varicose veins in pregnancy, prophylaxis of thrombosis and embolism.

CONTRAINDICATIONS:

Infections, eczema, skin ulcers of lower limbs, sever neuropathies, uncompensated heart failure, sever leg oedema, pulmonary oedema, pre-existing gangrene, hypodermitis, severe arterial circulatory disturbance.

AVAILABLE **STYLES**

VEGAL 1: 15-20 mmHg VEGAL 2: 22-28 mmHg

OPEN TOE AD AG AGH VEGAL 1 VEGAL 1 VEGAL 2 VEGAL 2 COLOUR

** AGG STYLE, LEFT OR RIGHT




COMPOSITION:

Cotton Polyamide Elastan

COLOUR







SIZES (cm)

S-M-L-XL-XXL

SIZE LEG CERCUMFERENCES (cm)

		S	М	L	XL	XXL	
	g	→ 48	→ 54	→ 60	→ 66	→72	
	d	→ 30	→ 35	→ 40	→ 45	→ 50	,
	b	17-20	20-23	23-26	26-29	29-32	