# SEQUENTIAL

## On the forefront of pneumatic compression is our 2000 series pumps. Designed to be user friendly, a digital

screen allows for easy operation. The pause button allows you to stop treatment and start again without skipping a beat; focus therapy allows for extra attention to a specific area; and the compliance meter keeps track of your hours of use. These gradient sequential pneumatic compression devices are used for the treatment of lymphedema, venous insufficiency, and various other circulatory disorders. Improve your quality of life starting today with the 2000 series pumps.





#### SC-2008-DL

#### Benefits Include:

- Manage pain
- Increase lymphatic and vascular flow
- Promote wound healing
- Reduce swelling
- Prevent infection
- Improve quality of life

#### **Indications**

- Lymphedema
- Lipedema
- Phlebolymphedema
- DVT prophylaxis
- · Chronic venous insufficiency
- Wound care

### SC-2004-DL

#### **Key Features:**

- 4-Chamber compression therapy
- Unilateral or bilateral capability
- Adjustable pressure between 10 mmHg and 120 mmHg
- Adjustable cycle time between 60 and 120 seconds
- Countdown timer or continuous mode
- Focus therapy
- Pause button
- Compliance meter
- 3 year warranty

Specifications:

Dimensions: 4.5" x 12" x 7.34"

11 cm x 30 cm x 18 cm

Weight: 3.3 lbs. 1,5 kgs

HCPCS: E0651

#### **Key Features:**

- 8-Chamber compression therapy
- Unilateral or bilateral capability
- Adjustable pressure between 10 mmHg and 120 mmHg
- Adjustable cycle time between 60 and 120 seconds
- Countdown timer or continuous mode
- Focus therapy
- Pause button
- Compliance meter
- 3 year warranty

Specifications:

Dimensions: 4.5" x 12" x 7.34"

11 cm x 30 cm x 18 cm

Weight: 3.65 lbs.

1,7 kg

HCPCS: E0651





See product manual for warnings, precautions, and contraindications.