

I ♥ my Legs

Gloria<sup>®</sup>  
1926



YASHODHAN ENTERPRISE

[www.connect2ye.com](http://www.connect2ye.com)



**MEDICAL COMPRESSION STOCKINGS IN SYNTHETIC, SOFT, LIGHTWEIGHT AND EASY TO WEAR**


Easyven® stockings, In soft and elastic synthetic yarn, combine fit and comfort without losing the therapeutic efficacy, Easy to put on and remove provide comfort and well-being throughout the day.

**INDICATION:**

moderate Chronic venous insufficiency (CVI) and after: healing of superficial ulcer, phlebitis syndrom, sclerotherapy treatment, ambulant selective varicose surgery, stripping, varices during pregnancy, prophylaxis of thrombosis/embolism.



Compression Class	Composition	Colour	Toe	Sizes
20-30 mmHg	66% PA 34% EA	Skin	Open	S - XL



**100% NATURAL COTTON FIBERS GENTLY CONTACT THE SKIN**



Stockings with high cotton content for a comfortable contact with the skin. Hypoallergenic & transpiratory, the stockings are well suited for those who suffer from dermatological problems. They are also ideal for warm & humid environments, for elderly people, & for people with allergies to synthetic fibers.

**INDICATION:**

Moderate symptoms of venous diseases such as aching heaviness & fatigue in legs, chronic venous insufficiency, post sclerotherapy, prophylaxis of thrombosis, moderate to severe varicosities, post phlebotic syndrome, lymphedema, stasis dermatitis due to venous insufficiency, post fracture or post surgical edema, edema from burns, varices during pregnancy, pregnant patients with previous phlebitis.

Compression Class	Composition	Colour	Toe	Sizes
23-32 mmHg	53% Cotton 30% EA 17% PA	Skin	Open	XS - XL



**GLORIA ULCERKIT® PRO MULTI-LAYER STOCKING SYSTEM**

Gloria UlcerKit Pro® is multi-layer stocking system to treat venous ulcers and is available in two different compressions: **UlcerKit 35 Strong (35-40 mmHg)**, **UlcerKit 30/20**

**INDICATION:**

**UlcerKit® Pro 35 (35-40mmHg)** : skin ulcers of the leg without venous insufficiency (eg post-traumatic, vasculitic), skin ulcers of the leg with mild to moderate venous insufficiency.

**CONTRAINDICATIONS:** Gloria UlcerKit Pro® should not be used in the presence of significant arterial disease with (ABI) of less than 0.8, Diabetic neuropathy.



Compression Class	Composition	Colour	Toe	Sizes
23-32 mmHg (1pcs) 18 mmHg (1 pcs) Thin understocking AD (1pcs)	78% PA 22% EA	Skin	Open	S - XL



**GLORIA MED® SLEEVE**




Elastic compression arm-sleeves GloriaMed® lightweight and comfortable, thanks to the use of fine synthetic yarn, provide excellent comfort to the superior limb.

**INDICATION:**

Arm-Sleeve are indicated in the treatment of primary and secondary lymphedema of the upper limbs.  
**Class2:** Prevention of edema after decongestion is achieved. Post-mastectomy fluid accumulation and edema following radiation therapy. Post-surgery and post-traumatic edema after basic treatment.

**CONTRAINDICATIONS:** Uncontrolled lymphedema, Inflammation, Open wounds, Infections, Arterial circulatory problem.

Compression Class	Composition	Colour	Mitten covered	Sizes
23-32 mmHg	67% PA 33% EA	Skin	Yes	S - XL





## SHEERNESS, ELEGANCE AND MEDICAL COMPRESSION

The use of state-of-the-art technologies and highly elastic yarns guarantee the best possible comfort and fit.

### INDICATION:

#### Class 1: (18-21 mmHg)

Prophylaxis of venous diseases, heaviness and fatigue in the legs, pain with or without varicosity, initial varices during pregnancy.



Compression Class	Composition	Colour	Toe	Sizes
18-21 mmHg	60% PA 40% EA	Skin	Open	XS - XL



AD

## A COMPLETE RANGE FOR DYNAMIC MEN

Preventive socks of high quality, elegance and comfort permit wearing them throughout the day whether at work or leisure. The cotton style provides increased transpiration and reduced allergies.

### INDICATION:

Designed specially to provide therapeutic benefits and compression in presence of heaviness and fatigue in the legs. Prophylaxis of venous disorders in patients at risk.

Compression Class	Composition	Colour	Toe	Sizes
18-20 mmHg	36% Cotton 14% EA 50% PA	Black, Blue	Closed	XS - XXL



## GRADUATED COMPRESSION SOCKS FOR MAN & WOMEN

Sitting or standing for long periods of time without moving your legs is a condition that cause a constriction of the veins and a slow return of venous blood to the heart. Swollen and painful legs are the first possible signs. The graduated compression socks for men and women "Travel & Work" contribute actively to the health of your legs.

### INDICATION:

#### 15-20 mmHg

Ideal for long trips by car, plane, train and for those who in everyday life maintains a static position for many hours, or with their legs crossed.

Compression Class	Composition	Colour	Toe	Sizes
15-20 mmHg	88% PA 12% EA	Black, Blue	Closed	S - XL



## COMFORT AND PROTECTION FOR DIABETIC AND SENSITIVE FEET

Diabesoft® is a line of socks indicated for people who suffer from diabetes, Athlete's foot, sensitive and delicate feet. The special cotton yarn with integrated silver has an antibacterial effect, anti-odor and regulate the temperature keeping the foot healthy, dry and protected. Without compression to not compress and avoid tourniquet effect. Without seams and with soft terry insole, to avoid skin abrasions and ensure a perfect fit and maximum comfort.

### INDICATION:

Recommended especially for people with diabetic sensitive and delicate feet, Athlete's foot.

Compression Class	Composition	Colour	Toe	Sizes
Without compression	98% Cotton (Silver Integrated) 2% EA	Black, Skin	Closed	S - XXL





**GLORIA MED**  
Thrombosis

## THROMBOEMBOLIC PROPHYLAXIS


A complete line of anti-embolism stocking. Comfortable and easy to wear thanks to a very soft fabric, well tolerate due to a normal transpiration of the skin. Compression of 18mmHg at the ankle.

### INDICATION:

Clinically designed to promote venous return for prevention of thrombosis before, during and after surgery and for medical conditions requiring prolonged bed rest. Before, during and after childbirth. Bedridden oncologic patients.

**CONTRAINDICATIONS:** Infections, eczema, Ulcers of lower limbs, severe neuropathies, uncompensated heart failure, severe leg edema, pulmonary edema, pre-existing gangrene, severe arterial circulation disorders.

Compression Class	Composition	Colour	Toe	Sizes
18 mmHg	65% PA 35% EA	White	Closed toe with inspection hole	S - XXL



Stockings are routinely tested by the Hohenstein Research Institute in Germany

## Extended Product Range



Conical fiber enables the use of guided laser energy to treat haemorrhoids of 3rd & 4th degree. This results in a shrinking & fibrotization without complications & anoderm survival with mucous membrane restoring normal anatomical structures. Specially designed conical tip of haemorrhoid fiber can be inserted directly in affected area without incision avoiding open wounds & excisions.

**Distal Tip:** Forward direction cone  
**Wave length:** 980nm to 1940nm  
**Fiber Core diameter:** 600µm  
**N.A.:** 0.20 up to 0.37  
**Length:** 2.5m



Radial fiber we can offer is the true 360° emitting probe serves for need of endo-venous laser treatment. The ink-free marking on the fiber sleeve eases positioning of the probe inside the vein. Radial fiber delivers a homogenous firing on to the vein wall without front-firing.

**Distal Tip:** 60° cone from fiber axis  
**Wave length:** 980nm to 1940nm  
**Fiber Core diameter:** 400µm & 600µm  
**N.A.:** 0.20 up to 0.49  
**Length:** 2.5 m



Bare Fiber we can offer a variety of fibers that are designed to meet different needs & applications. Our customers can choose between different core sizes as needed.

**Distal Tip:** Flat tip  
**Wave length:** 200nm to 2200nm  
**Fiber Core diameter:** 100µm up to 1000µm  
**N.A.:** 0.20 up to 0.49  
**Length:** 3 m



It is Unique Advanced Compression Bandage with Pocket and Velcro closure. The fine quality Advanced Compression bandage is manufactured as per international standards.

**Sizes:** 10 cm x 4 m, 15 cm x 4 m, 15 cm x 8 m

Notes:

.....

.....

.....

.....

Distributor Details:



**YASHODHAN ENTERPRISE**

**Sales & Admin Office :**

Office No. 8, Basement, Golden Heavens CHS, Near TJSB Bank, Kolbad Road, Khopat, Thane (W) 400601.

Contact : +91-22-25473636 / +91 70217 34147

Email : [info@connect2ye.com](mailto:info@connect2ye.com) [www.connect2ye.com](http://www.connect2ye.com)

### Graduated Compression Garments

Varicose Veins, DVT, Venous Insufficiency, Travellers, Diabetes, Sensitive Feet.

Medical Diode Laser Fibres (Bare, Radial, Conical)

Pre-powdered Latex examination Gloves

### Leasing Services

Diode laser machine for laser varicose vein, proctology, PLDD, Gynaecology.