



MEDICAL COMPRESSION STOCKINGS IN SYNTHETIC SOFT, LIGHTWEIGHT AND EASY TO WEAR

Easyven® stockings, in soft and elastic synthetic yarns, combine fit and comfort without losing the therapeutic efficacy. Easy to put on and remove provide comfort and well-being throughout the day.



COMPRESSION CLASSES:

EASYVEN 1: 15 - 20 mmHg **NEW**

EASYVEN 2: 20 - 30 mmHg

INDICATION:

EASYVEN 1: In presence of light venous disorders such as heaviness and fatigue in the legs, even during pregnancy. In non-serious post-traumatic situations, post-plaster.

EASYVEN 2: moderate chronic venous insufficiency (CVI) and after: healing of superficial ulcer, phlebitis syndrome, sclerotherapy treatment, ambulant selective varicose surgery, stripping. Varices during pregnancy, prophylaxis of thrombosis/embolism.

CONTRAINDICATIONS:


Infections, eczema, ulcers of lower limbs of arterial origins, severe neuropathies, uncompensated heart failure, severe leg edema, pulmonary edema, pre-existing gangrene, severe arterial circulation disorders.

AVAILABLE STYLES

EASYVEN

EASYVEN 1: 15 - 20 mmHg
EASYVEN 2: 20 - 30 mmHg

OPEN TOE



AD AGH AGG* AT

	EASYVEN 1	EASYVEN 1	EASYVEN 1	EASYVEN 1
	EASYVEN 2	EASYVEN 2	EASYVEN 2	EASYVEN 2
COLORS				

DETAILS

STAY-UP BAND:
plain with silicon drops

COMPOSITION:
Polyamide
Elastane

COLORS

34

Beige

* AGG: adjustable waist belt, left or right style

* Available from January 2022

SIZES (CM)

S - M - L - XL

NORMAL - MAX

SHORT - LONG

CIRCUMFERENCE MEASUREMENTS (CM)

	S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max
g	47 - 51	52 - 60	52 - 57	58 - 67	58 - 63	64 - 74	64 - 70	71 - 83
d	31 - 34	35 - 40	35 - 38	39 - 44	39 - 42	43 - 50	43 - 47	48 - 55
b	20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (CM)

	SHORT	LONG
<i>l</i> D	34 - 38	39 - 44
<i>l</i> G	62 - 71	72 - 82

